

# roto-mini

## PLUS

**ДИАМ**  
современная лаборатория

**www.dia-m.ru**  
заказ on-line

## FLEXIBLE VARIABLE SPEED ROTATOR FOR A RANGE OF APPLICATIONS

With speeds from 5 to 70 rpm, the Roto-Mini Plus variable speed rotator offers three operational mixing modes (Mix, Pause, and Rock) and accommodates standard tube sizes from 0.5 mL to 50 mL, making it a flexible addition to any lab.

### FLEXIBLE OPERATIONAL MODES

With speeds from 5 to 70 rpm, the Roto-Mini Plus variable speed rotator provides three different mixing modes to cover many applications, including:

- Mix – samples rotate with oscillations for increased mixing
- Pause – sample rotation includes periodic pauses
- Rock – samples gently rock in a see-saw manner.

### ADAPT TO MANY TUBE STYLES

The Roto-Mini Plus includes four holders that cover the following common tube sizes:

- 11 x 1.5/2.0 mL tubes
- 5 x 5/7 mL tubes
- 5 x 10/15 mL tubes
- 6 x 50 mL tubes/24 x 0.5 mL tubes



### ROTO-MINI PLUS

Description	Part Number
Roto-Mini Plus with tube holders, 115V, US Plug	36110100
Roto-Mini Plus with tube holders, 230V, EU Plug	36110110
Roto-Mini Plus with tube holders, 230V, UK plug	36110120

**000 «Диаэм»**

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

**www.dia-m.ru**

**С.-Петербург**  
+7 (812) 372-6040  
spb@dia-m.ru

**Новосибирск**  
+7 (383) 328-0048  
nsk@dia-m.ru

**Воронеж**  
+7 (473) 232-4412  
vm@dia-m.ru

**Екатеринбург**  
+7 (912) 658-7606  
ekb@dia-m.ru

**Йошкар-Ола**  
+7 (927) 880-3676  
nba@dia-m.ru

**Кемерово**  
+7 (923) 158-6753  
kemerovo@dia-m.ru

**Красноярск**  
+7 (923) 303-0152  
krsk@dia-m.ru

**Армения**  
+7 (094) 01-0173  
armenia@dia-m.ru

**Казань**  
+7 (843) 210-2080  
kazan@dia-m.ru

**Ростов-на-Дону**  
+7 (863) 303-5500  
rnd@dia-m.ru

